

Instruction Manual



QUICK GUIDE

- Wiring Diagram
Page 3-2
- Programming
the Display
Page 5-13
- Calibration
Chapter 6
- Troubleshooting
Chapter 7

Model 14 Dissolved Oxygen Transmitter

BROADLEY
TECHNOLOGIES

19 Thomas, Irvine, California 92618 USA
Call Toll-Free: 877.246.7900 (USA & Canada)
Phone: 949.452.1112 Fax: 949.452.1114
E-Mail: sales@broadleytech.com





Table of Contents

Chapter 1 : Introduction	1-1
1.1 Introduction.....	1-1
1.2 Principles of Operation.....	1-1
1.3 Main Characteristics	1-2
1.4 Technical Specifications.....	1-2
1.5 Dimensions.....	1-4
Chapter 2 : Description of the Analyzer	2-1
2.1 Transmitter	2-1
2.1.1 Presentation of the Transmitter.....	2-1
2.1.2 Transmitter Schematic.....	2-2
2.2 Application Fields	2-3
2.3 Standard Mounting Options	2-3
Chapter 3 : Installation of the Instrument	3-1
3.1 Unpacking.....	3-1
3.2 Inspection.....	3-1
3.3 Mounting	3-1
3.4 Location.....	3-1
3.5 Power Connection.....	3-1
3.6 Transmitter Start-Up.....	3-1
3.7 Adjusting the Display Contrast.....	3-1
3.8 Electrical Connections	3-2
3.9 Terminal Descriptions.....	3-3
3.10 Changing the Programming Language	3-4
3.11 Changing the Power Supply Frequency.....	3-5
Chapter 4 : Using the Instrument.....	4-1
4.1 Front Panel Keys	4-1
4.2 Displays 1 - 4.....	4-2
4.3 Description of the Function Keys.....	4-3
4.4 Icons.....	4-3
4.5 Enter or Modify a Value	4-4
Chapter 5 : Programming the Analyzer	5-1
5.1 Main Menu	5-1
5.1.1 CALIBRATION Menu	5-2
5.1.2 MAINTENANCE Menu.....	5-3
5.1.3 PROGRAMMING Menu.....	5-3
5.1.3.1 MEASURE Menu.....	5-4
5.1.3.2 ALARMS Menu	5-5
5.1.3.3 mA OUTPUTS Menu	5-8
5.1.3.4 RS485 Menu	5-10



Table of Contents (Continued)

Chapter 5 : Programming the Analyzer (Continued)	5-1
5.1.4 SERVICE Menu	5-11
5.1.4.1 AVERAGE Menu	5-12
5.1.4.2 DISPLAY Menu	5-13
5.1.4.3 CODE Menu	5-14
5.1.4.4 SOFT VERSION Menu	5-15
5.1.4.5 DEFAULT VAL. Menu	5-15
5.1.4.6 mA ADJUST Menu	5-16
5.1.4.7 CONFIGURATION Menu	5-16
Chapter 6 : Calibrating the Instrument	6-1
6.1 Calibrating the Temperature Sensor	6-1
6.1.1 Automatic Temperature Compensation	6-2
6.1.2 Manual Temperature Compensation	6-3
6.2 Calibrating the Measurement	6-4
6.2.1 Slope Calibration in Air with an Electric Zero	6-4
6.2.2 Slope Calibration in Air with Chemical Zero	6-5
6.2.3 Slope Calibration in the Process with an Electric Zero	6-6
6.2.4 Slope Calibration in the Process with Chemical Zero	6-7
Chapter 7 : Start-up and Troubleshooting	7-1
7.1 Start-up.....	7-1
7.1.1 Probe Connections	7-1
7.1.2 Main Power Supply Connection	7-1
7.1.3 Starting the Analyzer	7-1
7.2 Functional Troubleshooting	7-1
7.3 Troubleshooting the Electronics.....	7-3
Chapter 8: Error Messages	8-1
Appendix 1 : Pressure Conversion Table	A-1
Appendix 2 : Temperature Conversion Table	A-2
Appendix 3 : Default Values	A-3



Chapter 1: Introduction

1.1 Introduction


The Model 14 Dissolved Oxygen Analyzer is a single-channel analyzer for the measurement of dissolved oxygen in fermentation vessels and bioreactors, as well as in municipal applications.

1.2 Principle of Operation


The measurement of dissolved oxygen is based on the well-known Clark cell principle. An oxygen-permeable membrane isolates the electrodes from the sample, thus eliminating the need for sample conditioning. Other reducible or oxidizable ions do not interfere, because they cannot pass through the gas-permeable membrane.

A constant voltage supply powers two electrodes, maintaining each at a constant potential. A platinum work-

ing electrode (cathode) reduces the dissolved oxygen to hydroxyl ions :



A large silver counter electrode (anode) provides the oxidation reaction:



The reduction of oxygen is the current limiting reaction, thus making the cell current linearly proportional to the dissolved oxygen concentration.

Electrochemical reactions and diffusion rates are temperature-sensitive. The measuring cell, therefore, is equipped with a temperature sensor which allows an automatic temperature compensation.



1.3 Main Characteristics

- Range: 0 - 299.99% saturation; 0 - 29.99 ppm
- Calibration in the air
- Temperature compensation
- Programmable alarm levels, outputs on relays
- 4-20 mA, 0-20 mA analog outputs (standard)
- Wall-, panel- and pipe mounting

1.4 Technical Specifications

Electrical Characteristics

Power Supply	90 – 265 VAC , 50/60 Hz, Self-adjusting
Power Consumption	25 VA
Connections	Terminal blocks
Fuse	5x20 mm cartridge - T2AL - 250V

Analysis

Number of Inputs	1
Measuring Range	0 – 299.99% saturation
Accuracy	± 1% full scale
Reproducibility	± 0.1% of range
Response Time (90 %)	0 – 95% full scale, <45 seconds
Ambient Temperature	-20 – 60 °C (-4 ... 140°F)
Relative Humidity	10 – 90%



Transmitter

Display	Display in concentration units or % saturation Display of the direct cell current in nA Display of the sample temperature in °C/°F Programming via menus
Units	nA, ppm-mg/l, °C, °F, % saturation
Calibration	Electrical zero, chemical zero, slope calibration in the air, slope process calibration by comparison with a laboratory measurement
Analog Outputs	2 x 0/4 - 20 mA isolated from input signal, 800 ohms load maximum - Measurement or temperature - Mode : linear, bi-linear - Accuracy : 0.1 mA Alarms - Number : 4 - Functions : alarm - system alarm - timer - Hysteresis : 0 - 10% - Delay : 0 - 999 seconds - Breaking power : 250 VAC, 3A maximum 30 VDC, 0.5A maximum
Temperature Compensation RS485 (option)	Automatic in the range of 0 - 80 °C Speed : 300 - 9600 bauds Insulation : galvanic Station number : 32 maximum
Transmitter Protection Error Reports	IP 65 and NEMA 4X Cell current > 999 nA Sample temperature > 80 °C Ambient temperature > 80 °C Slope calibration error Zero calibration error (offset)

Electromagnetic Compatibility

Immunity Against Electromagnetic Interferences	EN 50082-2 and EN 50082-1
Electromagnetic Emission	EN 50081-1 and EN 50081-2
Low Voltage Standard	IEC61010-1

Materials

Working Electrode	Cathode : platinum
Counter Electrode	Anode : silver
Membrane Holder	PPS & silicone
Membrane	PTFE
Transmitter	Epoxy coated aluminum
Probe Body	Stainless steel 316L



Maintenance

Cleaning	Clean the instrument with a soft tissue. DO NOT use any aggressive agent.
----------	--



1.5 Dimensions

(Dimensions are in mm and inches)



Chapter 2: Description of the Analyzer

2.1 Transmitter

2.1.1 Presentation of the Transmitter



The electronic unit amplifies the signal of the amperometric measuring cell and converts it into a direct digital readout in ppm, mg/l or % saturation. The transmitter is comprised of the following items:


- Potentiostat which maintains the working electrode potential constant
- Amperometric measuring module
- Analog multiplexer
- Microprocessor unit

The analog multiplexer allows measurement inputs from the measuring cell, temperature sensor and internal checkpoints. Furthermore, the microprocessor operates the relays, the RS485 interface (optional) and the analog outputs.

The unit has an internal concentration auto-ranging feature and a microprocessor-operated calibration routine. The output of the potentiostat is monitored for possible overdriving of the potentiostat-output stage. This condition can occur when the connections to the measuring cell are due to an inoperable electrode or a defective reference electrode.



2.1.2 Transmitter Schematic:



1 : Programmable potentiostat	6 : Auxiliary input
2 : Polarization voltage amplifier	7 : Multiplexer
3 : Reversing switch for a 2 or 3 electrode operation.	8 : A/D converter
4 : Temperature measurement circuit	9 : Microprocessor
5 : Measurement circuit current amplifier	




2.2 Application Fields

Easy to use, install and program, this instrument is suitable for the following applications:

- Fermentation / biotech
- Municipal

2.3 Standard Mounting Options (using the red clamping bow)

The transmitter housing conforms to norm DIN 43700.




Panel mounting

Panel cut-out: 138 mm x 138 mm (5.4" x 5.4")

Front panel dimensions: 144 mm x 144 mm (5.7" x 5.7")

Included hardware: 2 flat head screws, 4 mm x 18 mm long, for panels up to 6 mm thick



Wall mounting

Included hardware: 2 pan head screws, 4 mm x 60 mm long for red clamping bow (a)


Additional required hardware: 2 flat head screws, 4 mm x 60 mm long (b)

Placed 3.15 " (80 mm) on center




Vertical pipe mounting

Included hardware: 2 pan head screws,
4 mm x 60 mm long for red clamping bow



Horizontal pipe mounting

Included hardware: 2 pan head screws,
4 mm x 60 mm long for red clamping bow



Chapter 3: Installation of the Instrument

3.1 Unpacking

The analyzer should be unpacked with great care. Watch for any loose accessories.

3.2 Inspection

The analyzer has been factory checked and tested prior to shipment. It is advisable, however, to inspect all parts immediately upon receipt for any damage which may have occurred during shipping. A damaged shipping container may indicate internal damage, which may not be immediately obvious. If there is any evidence of damage, keep the shipping container and refer to your local agent or to:

Broadley Technologies Corporation
19 Thomas, Irvine, CA 92618 USA
949-452-1112

Toll Free in the US and Canada 877-246-7900

3.3 Mounting

CAUTION!
Mounting should be done by qualified service personnel only. No power should be applied until the installation is complete.

3.4 Location

The analyzer should be located in an accessible site. The site should permit access for any inspection or maintenance operation.

3.5 Power Connection

For safety reasons, it is required to observe the precautions below:

1. The instrument should be connected to the power supply by means of a breaker located close to the instrument and clearly identified.
2. This breaker should switch off phase and neutral in case of electrical problems or to service the instrument. However, the earth ground must always be connected.

3.6 Starting the Transmitter

Before switching on the transmitter, make sure the site voltage corresponds to the instrument voltage indicated on the identification plate.

3.7 Adjusting the Display Contrast

If the display contrast is not sufficient, adjust the potentiometer P1 (blue color, see figure on page 3-2), which is located on the top left of the CPU board (after opening the enclosure).



Before servicing the instrument, confirm the power supply is "off".



3.8 Electrical Connections

■ Model 14

Do not switch the instrument on until completion of the installation.








An aluminum armor plate inside the Model 14 gives a detailed description of the different terminals and their connections:

⇒ *The Relays and Main Power Supply terminals represented on the right side are accessible by removing the armor plate.*






3.9 Terminal Descriptions

0/4 - 20 mA outputs galvanic insulation 	DESCRIPTION		CONNECTION
	0 - 20 mA or 4 - 20 mA (n°1) [+]		user
	0 - 20 mA or 4 - 20 mA (n°1) [-]		user
	0 - 20 mA or 4 - 20 mA (n°2) [-]		user
	0 - 20 mA or 4 - 20 mA (n°2) [+]		user
RS485 (option) 	RS485 Option		user
			user
amperometric module 	DESCRIPTION	COLOR	CONNECTION
	Temperature sensor [+]	black	temp +
	Temperature sensor [-]	blue	temp -
	Reference if using 3 electrodes	not used on the Model 14 oxygen analyzer	
	Anode	red	Anode
	Cathode	white	Cathode
	Auxiliary		N.C.
	External shield	braid	armor plate
	Internal shield	brown	GND
Preamplifier Supply 	Not used for the oxygen measurement		N.C.
Behind aluminum plate			
Mains 	Main power supply, 90...265 VAC 50/60 Hz or 24 V AC/DC (special version)		Ground
Amp module 	DESCRIPTION	CONNECTION	
	Alarm 1, simple contact	user	
	Alarm 2, simple contact	user	
	Alarm 3 or alarm system, simple contact	user	
	Alarm 4 or timer, simple contact	user	

⇒ Electrical connections should remain dry to ensure proper operation of the instrument. Check the creeping of the cables when opening the transmitter.
 ⇒ Shielded cables are recommended. This shielding should be connected to the earth central shielding.

3.10 Changing the Programming Language


The default programming language is English. To change the language, follow the procedure below (example for French):





3.11 Programming the Power Supply Frequency

The power supply frequency can be changed if necessary. This change occurs at the initial start-up and after resetting the instrument. Follow the procedure below.





Chapter 4: Using the Instrument

4.1 Front Panel Keys

The display may be programmed to indicate :

- Sample concentration
- Sample temperature
- Diffusion current
- Programming codes
- Programming features






Figure 4-1 : Front Panel


4.2 Displays 1 to 4 (live displays)




100% : Dissolved Oxygen measurement
23.2°C : Temperature measurement
O2 : Application
S1...S4 : Alarm status
(invisible if alarm is inactive)



Display of the parameters measured :
% Saturation
Temperature
Cell current



S1...S4 : alarm status
In this case relays S1 and S3 are active














Analog output allocation and level
Numeric and bargraph indication







4.3 Description of the Function Keys

The function keys below will be highlighted at the bottom of the screen :


	Modify a parameter
	Scrolling in a list of menus
	Return to the main display
	Display the main menu
	Display screen 2
	Display screen 3
	Display screen 4
	Validate the measure during calibration
	Confirm a command
	Decrease a value
	Increase a value

4.4 Icons

	Wait for instrument to reset
	Alarm system for relay S3
	Timer symbol : countdown for relay S4
	Controller symbol

4.5 Enter or Modify a Value

The highlighted digit can be modified with the key . Each digit can be validated by pressing ENTER. Repeat both operations for each digit.



WARNINGS!

Note 1: If you do not use the keyboard for at least 10 minutes, the instrument returns to the measuring mode.



Note 2: An access code may be required for the CALIBRATION, PROGRAMMING or SERVICE menus (see CODE menu).

It is possible to display a negative first digit "-"
It is possible to display a "." for the other digits.



Chapter 5: Programming the Transmitter

5.1 Main Menu




5.1.1 CALIBRATION Menus

- ⇒ Any calibration should follow the procedure below:
1. Configure the calibration characteristics in the “PROGRAMMING” menu.
 2. Perform the calibration via the “EXECUTION” menu.



An access code may be required if one has been programmed, see page 5-14 for CODE Menu.
The Temp. Calib. menu will not appear when instrument is set for manual temperature compensation.



Parameters		
DATE	mm/dd/yy	Date of the last calibration. The date programmed is not updated automatically.
S	x.xxx nA/ppm	Slope value
ΔT	x.x°C	Drift between the theoretical temperature (sensor curve) T_h and the temperature measured T_m : $DT = T_h - T_m$




5.1.2 MAINTENANCE Menu

<div style="background-color: black; color: white; text-align: center; padding: 2px;">MAINTENANCE</div> <p style="text-align: center; font-size: 1.2em;">21.0 ppm</p> <p>21.6°C 60.5 nA</p>	<p>Used for any maintenance operation in the instrument; the transmitter continues to display the variables measured.</p> <div style="background-color: black; color: white; padding: 5px;"> <p>The relay status is not modified. The analog output value depends on the configuration in the mA OUTPUTS/SPECIAL PROG. /MAINTENANCE menu.</p> </div>
--	--

5.1.3 PROGRAMMING Menu





An access code may be required.
See page 5-14 for CODE menu.





5.1.3.1 MEASURE Menu




Temperature Compensation		
TYPE	- Auto - Manual	Choice of a temperature measurement with automatic compensation or manual compensation  When manual temperature compensation is selected, the TEMP:CALIB. menu is no longer accessible!
TEMP.	- XX.X°C	In the case of manual compensation, enter the sample temperature



5.1.3.2 ALARMS Menu

Relays S1 through S4 may be allocated to the limit, alarm system or timer functions.



⇒ LIMIT FUNCTION:

The alarm relays are activated when the comparison between the measured value and the programmed limits meets the alarm function condition (up or down). The limits are programmed according to the following programming variables:

Limit Alarm Function			
<div style="border: 1px dashed black; padding: 5px;"> <p style="text-align: center;">ALARMS</p> <p>AFFECT.: Conc. LIM.: 0.001ppm DIR.: Down DELAY: 000s HYST.: 00% RELAY: NO</p> <p style="text-align: center;">Select </p> </div>	AFFECT	-Conc. - No -°C/°F	Programmed for a concentration or temperature limit or not active.
	LIM DIR.	XXXX -Up -Down	Enter a limit value Choice of the direction
	DELAY	XXxs	Time before the relay is executed (in seconds)
	HYST.	XX%	Definition of the hysteresis limit in % (10% max.). The hysteresis operates on only one side of the limit. The hysteresis is below the limit for the high alarm (up) and above the limit for the low alarm (down).
	RELAY	-NO -NC	Relay normally open or normally closed



⇒ SYSTEM ALARM FUNCTION:

Relays S3 and S4 may be used as a fault indicator. To control the faults traced by the analyzer, connect the specific relay to an external alarm system. The relay is activated as soon as a default appears.

In the case of a manual acknowledgment, the relay remains activated even if the default disappears. Press ENTER to deactivate the relay and the error message. In the case of an automatic acknowledgment, the relay and the error message are deactivated when the default disappears.

System Alarm			
	MODE	-No -Limit -System	The alarm S3 may be programmed to be a limit function (see paragraph above) or an alarm system function
	ACCEPT	-Auto -Manu.	In the case of an alarm system, choose between a manual (key ENTER) or an automatic acknowledgment
	RELAY	-NO -NC	Alarm S3 can be normally open or normally closed




⇒ TIMER FUNCTION:

Relays S3 and S4 may be programmed to a timer function.

Timer Function			
<div style="border: 1px dashed black; padding: 5px;"> <p style="text-align: center; background-color: black; color: white; margin: 0;">ALARM 4</p> <p>MODE: Timer</p> <p>INTERV: 1440mn</p> <p>IMPUL. Nb.: 5</p> <p>Ton: 005s</p> <p>Toff: 003s</p> <p>TmA: 05mn</p> <p style="margin-top: 10px;">Select ▲</p> </div>	MODE	-No -Limit -Timer	The Alarm 4 may be a limit (see parameters above) or a timer function
	INTERV	XXXXmn	Interval between 2 active cycles (in minutes)
	IMPUL.Nb:	X	Number of pulses during an active cycle
	Ton	XXXs	Adjustment of the relay active time (in seconds) for each pulse
	Toff	XXXs	Adjustment of the relay inactive time (in seconds) for each pulse
	TmA	XXmn	Hold time for the analog outputs after each cycle ⇒ <i>The analog output status depends on the configuration of the menu mA OUTPUTS/SPECIAL PROG./TIMER</i>

The measurement cycle lasts 4 seconds.


Example of a timer operating cycle :



5.1.3.3 mA OUTPUTS Menu


The analog output signals allow the transmission of the measurements from the analyzer to any external control system.

It is highly recommended to use shielded cables for the output signals. This shielding should be connected to the earth terminal on the armor plate.






Output 1/2		
AFFECT	- Conc. - nA - °C/°F	To set the choice of analog output to the concentration or temperature measurement
TYPE	- 0/20 - 4/20	Choice of the analog output type
MODE	- Lin - Dual	Choice between a linear or dual range (see drawing below)
LOWER	XXXX	Bottom scale value
MIDD.	XXXX	Mid-scale value (only in <i>dual mode</i>)
UPPER	XXXX	Top of the scale value
Special Prog.		
MODE	- last - preset - live	Characteristics of the analog output during calibration, alarm system, maintenance or timer active cycles. Display and output will be last stored value, a preset value, or a live measurement
VALUE	XX	Preset value (0 to 21 mA)
Test		Test the analog outputs in 1 mA increments (0 to 21 mA)





5.1.3.4 RS485 Menu (optional)


If the RS485 optional board is installed, program the parameters of the menu below. The optional RS485 board enables a connection between the analyzer and a digital communication system. The communication protocol is JBUS/MODBUS. Call Broadley Technologies for more information.



RS485		
N°	XX	Enter number (0 – 32)
BAUD	- 300 - 600 - 1200 - 2400 - 4800 - 9600	Transmission speed in baud
PARITY	- No - Odd - Even	Without parity bit With odd parity bit With even parity bit
BIT STOP	- 1 - 2	1 bit stop 2 bit stop



5.1.4 SERVICE Menu






5.1.4.1 AVERAGE Menu



An access code may be required.
See page 5-14 for CODE menu.

Program a moving average on the concentration measurement.


The measurement cycle lasts 4 seconds.



Average		
AVERAGE	X	Define the number of measurements to calculate the average (1-10)
Test		Display the difference between a measurement obtained with and without averaging



5.1.4.2 DISPLAY Menu




Display		
CONC	- % sat - ppm - mg/l	Choice of measurement units
TEMP.	- °C - °F	Choice of temperature units
PRESSURE	- mmHg - mbar - inHg	Choice of atmospheric pressure units
LANGUAGE	- F - GB - D - Sp - I	Choice of languages : - French - English - German - Spanish - Italian



5.1.4.3 CODE Menu

Protection codes may be programmed for access to the PROGRAMMING, CALIBRATION, and SERVICE menus. This code is deactivated by entering 0000.




Code		
CALIB.	XXXX	Access code to the "CALIBRATION " menu
PROG.	XXXX	Access code to the "PROGRAMMING" menu
SERVICE	XXXX	Access code to the "SERVICE" menu

To override the access code, press ESC and ENTER simultaneously to enter the menu selected.




5.1.4.4 SOFT ISSUE Menu

This menu displays the software version installed in the instrument.



5.1.4.5 DEFAULT Menu




Pressing YES will load the default values. The current programmed values, historic values and calibration parameters will be lost.

5.1.4.6 mA ADJUST Menu


The analog output signals are factory-adjusted (upper limit: 20mA). However, if the upper limit of one of the outputs, drifts, the span value can be adjusted with the mA ADJUST menu. Connect an ampere meter in series to the analog output terminals. Adjust the value until the ampere meter displays 20mA.

Note: The value displayed does not correspond to an actual mA value.



5.1.4.7 CONFIGURATION Menu

The operating frequency must be programmed to match the power supply. When the frequency is changed, the instrument is automatically reset.





Chapter 6: Calibrating the Instrument

NOTE

See Chapter 5 for programming commands.

REMARK

Any result (calibration or measurement) is always brought back to the reference temperature (25°C, 77°F). If the sample temperature is different from the reference temperature, manual or automatic temperature compensation must be used.


6.1. Calibrating the Temperature Sensor

The temperature sensor is located in the dissolved oxygen measurement probe. It is pre-set by the factory, but needs to be calibrated in the sample, on-site. The temperature sensor must be calibrated before calibrating the sensor's slope and zero.



6.1.1. Automatic Temperature Compensation


The sensor continuously measures the sample temperature. The concentration values are automatically calculated using the algorithms programmed in the transmitter.





6.1.2. Manual Temperature Compensation

Manual temperature compensation should be used only when the sample temperature is constant.




6.2. Measurement Calibration

The slope calibration can be done in air or water. The zero can be set with an electronic zero or a chemical zero. A chemical zero uses Nitrogen gas or oxygen-free solution.

6.2.1. Slope Calibration in Air with an Electronic Zero

Temperature compensation must be set to "auto."






6.2.2. Slope Calibration in Air with a Chemical Zero


1. Programming	
	<p>ZERO : Select a chemical zero Select slope in air</p> <p>SLOPE : Enter the atmospheric pressure in mmHg Enter the value of the ambient temperature</p>
2. Execution	
	<p>ZERO : Place the sensor in 99.999% Nitrogen gas.</p> <p>Press ENTER, the “CAL” message flashes and indicates the instrument is in calibration mode. Wait for the current to stabilize (approximately 1-2 min.) and press OK to validate the calibration. The oxygen concentration value flashes for 3 seconds. The instrument displays zero.</p>

6.2.3. Slope Calibration in Process with an Electronic Zero





6.2.4. Slope Calibration in Process with a Chemical Zero





Chapter 7: Start-up and Troubleshooting

7.1 Start-up

7.1.1 Probe Connection

See page 3-3 for the terminal connection.

7.1.2 Main Power Supply Connection

Remove aluminum plate inside the Model 14. Power should be connected to the “Mains” terminal indicated on the plate. See Page 3-3 for terminal connections.

7.1.3 Starting the Analyzer

When switched on, the analyzer performs an automatic test of its electronics and displays a first value. Wait until the measurement has stabilized (max. 30 min). Do not calibrate the analyzer before the temperature and the concentration are stable. (See Chapter 6.)

Note 1: Never shake the probe during a calibration, it may provoke an introduction of oxygen in the electrolyte and temporarily increase the cell current.

7.2 Functional Troubleshooting

PROBLEM: Excessive time for stabilization or no stabilization during air calibration.

Causes and solutions:

- A) There is an electrolyte leak (through the membrane). The current is too high because of an excessive penetration of oxygen.
⇒ **Change the membrane.**
- B) The electrolyte is contaminated.
⇒ **Change the electrolyte.**

PROBLEM: No significant current increase when the probe is in the air during calibration

Causes and solutions:

- A) The membrane has been torn.
⇒ **Change the membrane.**
- B) Cable or connections have been damaged.
⇒ **Check the connection to the transmitter. If the connection is correct check the connections of the probe connector.**

PROBLEM: Significant instability in measuring mode

Causes and solutions:

- A) There is humidity or water in the probe connector.
⇒ **Dry the connector. Check connections.**
- B) Incorrect connection.
⇒ **Check the connections of the transmitter to the probe.**
- C) Electromagnetic interferences close to the probe or transmitter cable.
⇒ **Find a better place for the cable and check the EMC level.**



PROBLEM: Inaccurate readings

Causes and solutions:

- A) The membrane permeability has changed (coating).
⇒ **Calibrate the analyzer and check if the concentration is back to normal.**
- B) Electrolyte contamination.
⇒ **Inspect probe components. Change the electrolyte.**
- C) Leak of electrolyte.
⇒ **Inspect probe components. Change the electrolyte.**
- D) Error during calibration or incorrect calibration.
⇒ **Calibrate again to check the parameters. If the error is confirmed, check the calibration current (too high, too low or unstable) and the concentration in the air. Refer to the problems described above.**
- E) The temperature calibration has not been performed.
⇒ **Check the temperature given by the transmitter and calibrate it (see page 6-2). If manual compensation has been selected, confirm the value is correct.**
- F) The sample temperature or pressure is out of the specifications.
⇒ **Change the probe location or modify the sample so that it meets the specifications.**

Other problems

Causes and solutions:

- A) The probe current is zero.
⇒ **There is no electrolyte in the probe (leak). There may be a short-circuit in the connection.**
- B) The probe current is negative.
⇒ **Connection problem to the anode circuit (loose contact).
⇒ Deposits on the anode.**
- C) The sample temperature corresponds to the specifications, but there is an error on the temperature.
⇒ **There may be a short-circuit on the temperature connection.**







7.3 Troubleshooting the Electronics

WARNING! Never attempt to service before disconnecting the instrument from the main power line.		
MALFUNCTION	POSSIBLE CAUSE	REMEDIES
No indication	No power; instrument is not connected correctly Defective fuse Instrument's power supply set for wrong line voltage Ribbon cable connecting power with CPU board not properly plugged in Connection between CPU board and measurement module loose Short circuit in power supply board Hardware is defective	Check for power, then check if connected properly Check fuse Check jumpers on power supply board for correct voltage settings Verify that the ribbon plugs are connected properly Check plug connections Visually check power supply board for shorts Call the Service Technician
LCD displays undefined characters	Malfunctioning CPU board or processor CPU hardware	Using the Instruction Manual, program the instrument to load the default values RESET the instrument by temporarily interrupting the power (5-10 seconds) Call the Service Technician
Keyboard does not operate; all keys are inactive	CPU malfunctioning, external interferences	If there is no response, RESET the instrument by temporarily interrupting the power (5-10 seconds). Check each key again. If there is no change, call the Service Technician.
Measurement is not correct	Instrument was programmed incorrectly System, including probe, not calibrated correctly Probe connected wrong Probe malfunctioning, possibly incompatible with the application CPU board is defective	Recheck programmed parameters. Do they agree with the probe's characteristics? Calibrate the whole system (with probe connected) Recheck all probe connections Visually check the condition of the probe. Is the application within the probe's specifications? If error persists, call the Service Technician
Measurement is not stable	Faulty probe Probe connected wrong	Check the condition of the probe. Is it contaminated? Verify the probe is connected properly



MALFUNCTION	POSSIBLE CAUSE	REMEDIES
	Interferences Cable shield is not connected Defective CPU board	Are there any sources of potential interference, chemical, external, temperature, pressure, etc.? Check and connect If problem persists, call the Service Technician
Temperature measurement is not correct	Probe connected wrong Temperature not calibrated CPU board is defective	Check probe connection Calibrate for temperature. Also check the T.C. for correct resistance (ohmmeter) If problem persists, call the Service Technician
Display reading static; cannot be changed in any way	CPU board is malfunctioning and/or some other transmitter hardware is defective	Check probe connection Initiate a software RESET If Steps 1-3 do not remedy the problem, make a cold start (RESET): interrupt the power for 5-10 seconds Reprogram the instrument If problem persists, call the Service Technician
Relays not energized	Instrument was programmed incorrectly Hardware is defective	Check whether the correct relay parameters and set points have been programmed Check that the programmed set points are compatible with the programmed measuring range Check the relay characteristics for proper functioning using an ohmmeter If problem persists, contact the Service Technician
Wrong output current, output current remains locked at 0 or 20 mA	Instrument was incorrectly programmed Connection of the Model 14 with peripherals (recorder, etc.) are faulty, loose or defective Hardware is defective	Check the programmed output-current parameters Check the cables Compare the measured value with the output-current range If problem persists, call the Service Technician
Polarization voltage incorrect	Wrong configuration (3 electrode mode)	Check that the switches under the amperometric board are correctly positioned ("on")
Wrong programming	Check programmed parameters	

Chapter 8: Error Messages

During failure, the measurement is replaced by dashes “- - -.”

Error Message	Description
Error messages during a measurement	
	<p>The concentration value exceeds the limits. Check the current value and the calibration parameters.</p> <p>Electrode protection If the measuring range displays “concentration too high” for more than 2 minutes, the error message “STAND-BY” appears. From this moment on, the electrode polarization voltage is switched off for 10 minutes before the voltage is switched on again.</p> <p>In STAND-BY the user may switch the voltage on at any time by pressing ENTER.</p>
	<p>The concentration value exceeds the limits. Check the current value and the calibration parameters.</p>
	<p>The sample temperature exceeds the limits. Check the cable polarity. Check for a short-circuit or an open circuit.</p>
	<p>The current value exceeds the limits. Check the electrode (electrolyte and membrane).</p>



	<p>The current value exceeds the limits. Confirm there is no short-circuit in the measuring line. Check the polarization voltage.</p>
<p>Error messages during calibration</p>	
	<p>The temperature difference between the calibration temperature and the sensor's theoretical exceeds the limit programmed. Limits : $\pm 20^{\circ}\text{C}$</p>

Appendix 1: Pressure Conversion Table

hPa or mbars	mm Hg	mm H ₂ O	psi	inches Hg	inches H ₂ O
920.00	690	9503.60	13.34	27.16	374.16
925.00	693.8	9555.25	13.42	27.31	376.19
930.00	697.5	9606.90	13.49	27.46	378.22
935.00	701.3	9658.55	13.56	27.61	380.26
940.00	705	9710.20	13.63	27.76	382.29
945.00	708.8	9761.85	13.71	27.90	384.32
950.00	712.5	9813.50	13.78	28.05	386.36
955.00	716.3	9865.15	13.85	28.20	388.39
960.00	720	9916.80	13.92	28.35	390.43
965.00	723.8	9968.45	14.00	28.50	392.46
970.00	727.5	10020.10	14.07	28.64	394.49
975.00	731.3	10071.75	14.14	28.79	396.53
980.00	735	10123.40	14.21	28.93	398.56
985.00	738.8	10175.05	14.29	29.08	400.59
990.00	742.5	10226.70	14.36	29.23	402.63
995.00	746.3	10278.35	14.43	29.38	404.66
1000.00	750	10330.00	14.50	29.53	406.69
1005.00	753.8	10381.65	14.58	29.68	408.73
1010.00	757.5	10433.30	14.65	29.82	410.76
1015.00	761.3	10484.95	14.72	29.97	412.79
1020.00	765	10536.60	14.79	30.12	414.83
1025.00	768.8	10588.25	14.87	30.27	416.86
1030.00	772.5	10639.90	14.94	30.42	418.89

Appendix 2: Temperature Conversion Table

- Conversion from °C into °F : $^{\circ}\text{F} = 1.8 * ^{\circ}\text{C} + 32$
- Conversion from °C into °K : $^{\circ}\text{K} = ^{\circ}\text{C} + 273.15$

°C	°F	°K
0	32	273.15
1	33.8	274.15
2	35.6	275.15
3	37.4	276.15
4	39.2	277.15
5	41	278.15
6	42.8	279.15
7	44.6	280.15
8	46.4	281.15
9	48.2	282.15
10	50	283.15
11	51.8	284.15
12	53.6	285.15
13	55.4	286.15
14	57.2	287.15
15	59	288.15
16	60.8	289.15
17	62.6	290.15
18	64.4	291.15
19	66.2	292.15
20	68	293.15
21	69.8	294.15
22	71.6	295.15

°C	°F	°K
23	73.4	296.15
24	75.2	297.15
25	77	298.15
26	78.8	299.15
27	80.6	300.15
28	82.4	301.15
29	84.2	302.15
30	86	303.15
31	87.8	304.15
32	89.6	305.15
33	91.4	306.15
34	93.2	307.15
35	95	308.15
36	96.8	309.15
37	98.6	310.15
38	100.4	311.15
39	102.2	312.15
40	104	313.15
41	105.8	314.15
42	107.6	315.15
43	109.4	316.15
44	111.2	317.15
45	113	318.15

Appendix 3: Default Values

Calibration

CONC. CALIB.
ZERO CALIB.
 Type : Elec. Auto
SLOPE
 Type : Air
 PO: 0760.0 mm Hg

PARAMETERS
 DATE :01/01/98
 ZERO: 0.00 nA
 S: 3.000 nA/ppm
 ΔT: 0.0°C

Programming

MEASURE

TEMP. COMP.
 TYPE : Auto

ALARMS

ALARMS S1/S2/S4
 AFFECT. : Conc.
 DIR. : Down
 LIMIT: 0.000 ppb
 DELAY : 000 s
 HYST. : 00%
 RELAY : NO

ALARM S3
 AFFECT. : System
 ACQUIT : Auto
 RELAY : NC

mA OUTPUTS

OUTPUT 1
 AFFECT. : Conc.
 TYPE : 4-20
 MODE : Lin.
 LOW : 0.000 ppb
 UP : 01.00 ppm

OUTPUT 2
 AFFECT. : Conc.
 TYPE : 4-20
 MODE : Lin.
 LOW : 0.000 ppb
 UP : 01.00 ppm

SPECIAL PROG.

MAINTENANCE
 MODE : Last

CALIBRATION
 MODE : Last

SYST. ALARM
 MODE : Last

TIMER
 MODE : Last

RS485

N^o : 00
 BAUD : 9600
 PARITY : No
 STOP BIT : 1



Service

AVERAGE

AVERAGE : 1

DISPLAY

DISPLAY

CONC. : ppb/ppm

TEMP. : °C

PRESSURE : mmHg

LANGUAGE : GB

CODE

CODE

CALIB. : 0000

PROG. : 0000

SERVICE : 0000

CONFIGURATION

CONFIGURATION

FREQ. : 60 Hz



Broadley Technologies Corporation

19 Thomas
Irvine, California 92618
USA

Corporate Headquarters:
Toll-Free: 877.246.7900
Tel: 949.452.1112 Fax: 949.452.1114
E-Mail: sales@biotechcatalog.com

Copyright © 2000 Broadley Technologies Corporation
All Rights Reserved