

New Single-Use Ultrasonic Flowmeter

PendoTECH, USA have recently added a range of four LEVIFLOW[®] single-use ultrasonic flowmeters to their sensor product line. The flowmeter uses ultrasonic transit time, flow techniques to provide an obstruction free flow path, with no moving parts. The flowmeters find many different applications in the biopharm industry, including TFF filtrate flow, Normal flow/depth filtration and Chromatography.

Figure 1 illustrates the operating principle of the ultrasonic flowmeter.



operating principle of the ultrasonic flowmeter.

Two piezoelectric transducers, mounted in the sensor housing, generate and receive an ultrasonic wave. The wave going in direction of the flow (with-stream wave) is accelerated and the wave going against the flow direction (against-stream wave) is slowed down. The two waves are processed by the PendoTECH Leviflow[®] Sensor Monitor. The difference of the transit time of both waves is proportional to the velocity of the fluid. The monitor has a digital LED display for the flow reading. It also has both a 4-20mA analogue output and a digital frequency output. These outputs facilitate interface of the monitor to other systems for process control and data acquisition.

The standard configuration of the LEVIFLOW single-use flowmeters (Figure 2) consists of a flow sensor and the PendoTECH LEVIFLOW Sensor Monitor with a digital signal processor (DSP) for processing the sensor signals.



Single Use Ultrasonic Flowsensor

SYSTEM BENEFITS

- High precision (1%) flow measurement
- Flow sensor can be used as single use to prevent cross-contamination or re-used.
- No moving parts -> no particle generation
- Improved bubble robustness due to DSP (digital signal processing) technology
- Easy integration into OEM equipment
- FDA, USP-VI, BSE/TSE and Animal free wetted materials
- Gamma sterilizable polypropylene (PP).

The contents of this publication are presented for information purposes only, and while every effort has been made to ensure their accuracy, they are not to be construed as warranties or guarantees, express or implied, regarding products or services described herein or in their use or applicability. All sales are governed by our terms and conditions, which are available on request. We reserve the right to modify or improve the design or specification of such products at any time.

© 2024 Broadley-James Corporation. All rights reserved. Visit www.broadley-james.com/trademarks for trademark information.

TMP-BF-102101

New Single-Use Ultrasonic Flowmeter

The sensors are available in 4 different flow ranges to cover the majority of biopharm applications in process development and pilot plant. Full scale of 0.8 lpm, 8.0 lpm, 20 lpm and 80 lpm, with a typical accuracy of +/- 1% (see data sheet for full details).

Download [Single Use Flowmeter data sheet](#)

Broadley-James Ltd is the exclusive distributor for PendoTECH sensors in the UK and for Automation Systems throughout Europe. If you have a measurement or control application in upstream or downstream bioprocessing contact us now on sales@broadleyjames.co.uk

The contents of this publication are presented for information purposes only, and while every effort has been made to ensure their accuracy, they are not to be construed as warranties or guarantees, express or implied, regarding products or services described herein or in their use or applicability. All sales are governed by our terms and conditions, which are available on request. We reserve the right to modify or improve the design or specification of such products at any time.

© 2024 Broadley-James Corporation. All rights reserved. Visit www.broadley-james.com/trademarks for trademark information.

TMP-BF-102101



North America and Pacific

Email: helpdesk@broadleyjames.com
Web: www.broadleyjames.com

19 Thomas, Irvine CA, 92618 USA

Phone: 949-829-5555
Tollfree: 800-288-2833
Fax: 949-829-5560

United Kingdom and EU Countries

Email: sales@broadleyjames.co.uk Phone: +44 (0)1525 862518
Web: www.broadleyjames.eu Fax: +44 (0)1525 862811

Wrest Park, Silsoe Beds MK45 4HS, UK