

Description

The SingleSense®, USP Class VI pH Electrode is gamma sterilizable with standard European connectors and dimensions.

Purchase Information

Part Number	SU200-225-S8
Description	pH Electrode, S8, 225 mm, SingleSense®

Classification

Model Type	SU200
Product Type	pH Electrode

Specifications

Measurement

Type	pH
Principle	Potentiometric
Range	2–11 pH
Drift	< 0.05 pH per week

Operational Range

Temperature	0°–80° C [32°–176° F]
Pressure	0–2327 mmHg [0–45 psig]
Flow Rate	Not applicable

Signal

Input	Not applicable
Output	Analog: ± 414 mV

Sterilization

Autoclavable	No
Steam Sterilizable	No
Gamma Sterilizable	Yes, up to 50 kGy

Connectors and Cables


Connector	S8
Cable Length	Not applicable

Process Interface

Insertion Method	Pg13.5 Threaded Retainer Ring
Insertion Length	225 mm [8.85 inches]

Wetted Materials

Lead-free Glass; Alumina Silicate



SingleSense® pH Electrode
(12 x 225 mm)

Packaging and Regulations

Packaging Information

Package Type Cardboard box

Package Dimensions 40.64 x 10.16 x 5.72 cm [16 in x 4 in x 2.25 in]

Package Weight < 0.5 kg [1 lb]

Double Bagged No

Cleanroom Bagged No

Cleanroom Class..... Not applicable

Storage Conditions..... Temperature: 5°–50° C [41°–122° F]
Humidity: 5–95 RH

Included Accessories Not applicable

Quality Certificate..... Included

Shelf Life 30 months

Shelf Life Post Gamma Sterilization 30 months

Origination and Tariff Code

Harmonized Tariff Trade Code 9027803100

Country of Origin USA

MFR Location California, USA

Purchasing Information

Unit of Sale..... Each

Minimum Purchase Quantity..... Each

Regulatory and Compliance

CE Compliant Not applicable

RoHS Compliant Yes

REACH Compliant Not applicable

CFR 21-11 Compliant Not applicable

BSE/TSE/ADI Compliant..... Yes

USP Class VI Compliant Yes, all wetted parts

FDA Material Compliant..... Yes, all wetted parts

The contents of this publication are presented for information purposes only, and while every effort has been made to ensure their accuracy, they are not to be construed as warranties or guarantees, express or implied, regarding products or services described herein or in their use or applicability. All sales are governed by our terms and conditions, which are available on request. We reserve the right to modify or improve the design or specification of such products at any time.



North America and Pacific

Email: sales@broadleyjames.com
Web: www.broadleyjames.com

Phone: 949-829-5555
Toll-free: 800-288-2833
Fax: 949-829-5560

19 Thomas, Irvine, CA, 92618 USA

United Kingdom and EU Countries

Email: sales@broadleyjames.co.uk
Web: www.broadleyjames.eu

Phone: +44 (0)1525 862518
Fax: +44 (0)1525 862811

Wrest Park, Silsoe, Bedford MK45 4HS, UK