

What are the benefits of using a single-use sensor?

A single-use sensor is typically shipped gamma sterilized and ready to use, eliminating the costly and time-consuming cleaning, maintenance, and operator installation required for reusable sensors. Single-use sensors reduce downtime between batches while providing process measurement and control as good or better than traditional autoclavable and sterilizable sensors.

The contents of this publication are presented for information purposes only, and while every effort has been made to ensure their accuracy, they are not to be construed as warranties or guarantees, express or implied, regarding products or services described herein or in their use or applicability. All sales are governed by our terms and conditions, which are available on request. We reserve the right to modify or improve the design or specification of such products at any time.

© 2026 Broadley-James Corporation. All rights reserved. Visit www.broadley-james.com/trademarks for trademark information.

FAQ-5363
TMP-BF-102101