

What is downstream processing?

Downstream processing typically refers to the post-harvest processing of the clarified production media. These downstream process steps include capture/recovery of the desired biologic, viral inactivation, concentration and purification, final formulation, and final filling of vials for packaging and distribution.

The contents of this publication are presented for information purposes only, and while every effort has been made to ensure their accuracy, they are not to be construed as warranties or guarantees, express or implied, regarding products or services described herein or in their use or applicability. All sales are governed by our terms and conditions, which are available on request. We reserve the right to modify or improve the design or specification of such products at any time.

© 2026 Broadley-James Corporation. All rights reserved. Visit www.broadley-james.com/trademarks for trademark information.

FAQ-5365
TMP-BF-102101