

Why is my FermProbe turning brown?

The longer the electrode is exposed to high temperatures (above 100 degrees C) the darker it becomes. Likewise, the higher the temperature, the darker the color. This is normal. If the electrode can be calibrated, it's OK to use it.



The contents of this publication are presented for information purposes only, and while every effort has been made to ensure their accuracy, they are not to be construed as warranties or guarantees, express or implied, regarding products or services described herein or in their use or applicability. All sales are governed by our terms and conditions, which are available on request. We reserve the right to modify or improve the design or specification of such products at any time.

© 2022 Broadley-James Corporation. All rights reserved. Visit www.broadley-james.com/trademarks for trademark information.

TMP-BF-102101